



Recombinant Uncharacterized protein K06H7.7 (K06H7.7)

Product Code	CSB-MP339465CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34515
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSFSAGESSA ATGDPQDFNK FFGEINSPSE AAKTIKAIK PSVALPESS CASYCASWVN DQSFEILPTD EENRQTEFSS SPLSDSVHSC LGFNRKKERK RRTSGDIEQT DKRKRTEERP SDKITFASNV FDNLARVTAE RDNPEKLSKY KNIPGVGTRE NIDEGISDIN QPSSKSSGVK SKYVEVMSLF QKSIQKFSIM NLDYMYTFDW NHYSKKNIQY VVFTKVCQSH APGHFVEMLK EDCSCRDGCS ADLAVKVRLT VPLDLSDIPI PCSEDIQIHV SFLADDLLER HISSNVNRSS NRDDLDFSR R TCEKMYDLYQ KMMRGQSSRM YTEWFDKCKI FVKNFLFSKS LHLYHGNVIR RKIEESQTLN KRRCFIFAQI FDLE
Source	Mammalian cell
Target Names	K06H7.7
Protein Names	Recommended name: Uncharacterized protein K06H7.7
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.