



Recombinant Putative zinc finger protein R05D3.3 (R05D3.3)

Product Code	CSB-BP339469CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34536
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MESPNMRKRL IEMESSLEVK KTATTVKQEP LETSSPLLAS MISTAPNPPP TSPHALNYLT LLAAFRQFLH GYIQAGGDEN LANTLGNYVA ALPMPPNSLN PIGPLNAPKS PFKQQQKFHS IIKAVKREPV SEPAQNAGSG DSSVTLSSAL GEDAVPQNRG SCTSMKPVPR FVMLPVVAAQ SASPSASGPT SEIQQLREAA FGRYKNVLCV ICNEWICSRN RKNHIEAHLN YRPYKCSACS YARRREIFVD QHIRTQHKG V EGVVMLSVDV LHVAMEVDRL AECEVTRTRK IIDTMQEKKD GDFGENKDFD EKALQMMLAE EAENKVVVIE SVSQVRPKVA NYHRRQR TKV LKKIYD TDVA KQTELKIVKD EEGGVPTMSI KLEEGFEINL EDLMKMVGAG SSITDSNEPG PSEIKKELAE V
Source	Baculovirus
Target Names	R05D3.3
Protein Names	Recommended name: Putative zinc finger protein R05D3.3
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.