



# Recombinant Uncharacterized protein R10E11.5 (R10E11.5)

<b>Product Code</b>	CSB-BP339470CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34549
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASSAGRDKL RSRGQRVFAF GSSTPRDL SH MSKVPQNL RN YDAKSVDKSP SDAPSLHDYI VKARSL SREA CDSRNQFVFG SSTPRTL AHL DKIPHKQRVY DAKIPKKNTT HSDFKAAPIR FNAPPVPITK PAKKEENRHV ITSKDDDIAS EPDIVQDREE FMNEMKKHKI ELEKKKSLGK NDSKTTIEIP TPKAAETQQE HVKFAHAPEI AGNSQKAVRI QSETQEEAPV AVKNLSANKM NDQIEVSQLM NEITPESVPA VESLDNQKNA NVVGDLLAKV QKVADDTIDK SKTTVAADVA KMSGALQKAE EEVVQTIDQT VKNIKSNVNE VKKDVEKNIA EKVDDITKEL EKSAKSLEET TDKIGSKIDN TSQAIKSNLE EASLKTEKSV NDAVKSGEKL IGDLASEAKK TFDSAEQKLD AKFSKIGSTA DSTLEDVKNG LRHF
<b>Source</b>	Baculovirus
<b>Target Names</b>	R10E11.5
<b>Protein Names</b>	Recommended name: Uncharacterized protein R10E11.5
<b>Expression Region</b>	1-444
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.