



Recombinant Uncharacterized protein R10E11.5 (R10E11.5)

Product Code	CSB-EP339470CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34549
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MASSAGRDKL RSRGQRVFAF GSSTPRDL SH MSKV PQNLRN YDAKSVDKSP SDAPSLHDYI VKARSL SREA CDSRNQFVFG SSTPRTL A HL DKIPHKQRVY DAKIPKNTT HSDFKAAPIR FNAPPVPITK PAKKEENRHV ITSKDDDIAS EPDIVQDREE FMNEMKKHKI ELEKKKSLGK NDSKTTIEIP TPKAAETQQE HVKFAHAPEI AGNSQKAVRI QSETQEEAPV AVKNLSANKM NDQIEVSQLM NEITPESVPA VESLDNQKNA NVVGDLLAKV QKVADDTIDK SKTTVAADVA KMSGALQKAE EEVVQTIDQT VKNIKSNVNE VKKDVEKNIA EKVDDITKEL EKSAKSLEET TDKIGSKIDN TSQAIKSNLE EASLKTEKSV NDAVKSGEKL IGDLASEAKK TFDSAEQKLD AKFSKIGSTA DSTLEDVKNG LRHF
Source	E.coli
Target Names	R10E11.5
Protein Names	Recommended name: Uncharacterized protein R10E11.5
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.