



# Recombinant Arabidopsis thaliana Heterodimeric geranylgeranyl pyrophosphate synthase large subunit 1, chloroplastic (GGPPS1)

<b>Product Code</b>	CSB-YP339503DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34802
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSVV TKEDNLRQSE PSSFDFMSYI ITKAELVNKA LDSAVPLREP LKIHEAMRYS LLAGGKRVRP VLCIAACELV GGEESTAMPA ACAVEMIHTM SLIHDDLPCM DNDDLRRGKP TNHKVFGEDV AVLAGDALLS FAFEHLASAT SSDVVSPVRV VRAVGELAKA IGTEGLVAGQ VVDISSEGLD LNDVGLEHLE FIHLHKTAAL LEASAVLGAI VGGGSDDEIE RLRKFARCIG LLFQVDDIL DVTKSSKELG KTAGKDLIAD KLTYPKIMGL EKSREFAEKL NREARDQLLG FDSDKVAPLL ALANYIAYRQ N
<b>Source</b>	Yeast
<b>Target Names</b>	GGPPS1
<b>Protein Names</b>	Recommended name: Heterodimeric geranylgeranyl pyrophosphate synthase large subunit 1, chloroplastic Short name= GGPP synthase 1 Short name= GGPS1 EC= 2.5.1.- Alternative name(s): (2E,6E)-farnesyl diphosphate synthase 1 Di
<b>Expression Region</b>	57-371
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.