



Recombinant Arabidopsis thaliana Floral homeotic protein APETALA 1 (AP1)

Product Code	CSB-EP339569DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35631
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGRGRVQLKR IENKINRQVT FSKRRAGLLK KAHEISVLCD AEVALVVFSSH KGKLFESTD SCMEKILERY ERYSYAERQL IAPESDVNTN WSMEYNRLKA KIELLERNQR HYLGEDLQAM SPKELQNLEQ QLDTALKHIR TRKNQLMYES INELQKKEKA IQEQNSMLSK QIKEREKILR AQQEQWDQQN QGHNMPPPLP PQQHQIQHPY MLSHQSPFL NMGGLYQEDD PMAMRRNDLE LTLEPVYNCN LGCFAA
Source	E.coli
Target Names	AP1
Protein Names	Recommended name: Floral homeotic protein APETALA 1 Alternative name(s): Agamous-like MADS-box protein AGL7
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.