



Recombinant *Saccharomyces cerevisiae* Protein PRY2 (PRY2)

Product Code	CSB-EP339637SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36110
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	AP VAVTVTQHVH QAATVVVQGI VRVENGQTLT TFITKGTQTA SASPVATTSA PIVVANAQVD SIATSVIQES AVVAESATFE ESSTETSEAF STATATIQAV QTSASATQDD VTTTLTSSSTQ PTSTTTPTTT TTSPTTTTSP TTTASPTTTA SPTTATTTQS TASSTQSSSS DFSTSMVNEH NTKRALHKDT GSLTWSDTLA TYAQNYADSY DCSGNLVHSG GPYGENLALG YGTTGSVDAW YNEITSYDYS NPGFSESAGH FTQVVWKGTS EVGCGLKSCG GEWGDYIICS YKAAGNVIGE FADNV MPLA
Source	E.coli
Target Names	PRY2
Protein Names	Recommended name: Protein PRY2 Alternative name(s): Pathogen related in Sc 2
Expression Region	19-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.