



Recombinant *Saccharomyces cerevisiae* DASH complex subunit SPC34 (SPC34)

Product Code	CSB-EP339641SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36131
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGESLDRCID DINRAVDSMS TLYFKPPGIF HNAILQGASN KASIRKDITR LIKDCNHDEA YLLFKVNPEK QSVSRRDGKE GVFDYVIKRD TDMKRNRRLG RPGEKPIIHV PKEVYLNKDR LDLNNKRRRT ATTSGGGLNG FIFDTDLIGS SVISNSSSGT FKALSAVFKD DPQIQRLLYA LENGSVLMEE ESNNQRRKTI FVEDFPTDLI LKVMAEVTDL WPLTEFKQDY DQLYHNYEQL SSKLRFIKKE VLLQDDRLKT MSQYHPSSSH DVAKIIRKEK DEIRRLEMEI ANLQE
Source	E.coli
Target Names	SPC34
Protein Names	Recommended name: DASH complex subunit SPC34 Alternative name(s): 34 kDa spindle pole component protein Outer kinetochore protein SPC34
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.