



# Recombinant Schizosaccharomyces pombe Ribonucleoside-diphosphate reductase small chain (suc22)

<b>Product Code</b>	CSB-MP339706SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P36603
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGLEHLEEFY YPKEHGEEVE YDSEQGVRKI YVKSİKETFN FDNVSEEEKQ EGGDYYLGKK EDELDEVVLR PNPFRFVLFY IKYHEIWQFY KKAASFVTA EEIDLSKDLV DWDNKLNADE RYFISTVLAY FAASDGIVNE NLLERFSSEV QIPEARCVYG FQIMIENIHS ETYSLLLDY IREPKEKQRH FDAILTMGSI KAKAKWALRW INDEDSTYAI RLVAFAAVEG IFFSGSFASI FWLKKRGLMP GLTFSNELIC RDEGLHTDFA CLMFSLKHR PGRKVVEAII VEAVDIEKEY FTDALPVSLG MNKDLMCQY IEFVADRLLV ALGNDKYYNV TNPFDPMENI SLAGKTNFFE KKVSDYQIAG VMSGTKRAEK DDHTFTIDED F
<b>Source</b>	Mammalian cell
<b>Target Names</b>	suc22
<b>Protein Names</b>	Recommended name: Ribonucleoside-diphosphate reductase small chain EC=1.17.4.1 Alternative name(s): Ribonucleotide reductase small subunit
<b>Expression Region</b>	1-391
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.