



# Recombinant Escherichia coli Blue copper oxidase CueO (cueO)

<b>Product Code</b>	CSB-EP339713ENV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P36649
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AE RPTLPIPDLL TTDARNRIQL TIGAGQSTFG GKTATTWGYN GNLLGPAVKL QRGKAVTVDI YNQLTEETTL HWHGLEVPGE VDGGPQGIIP PGGKRSVTLN VDQPAATCWF HPHQHGKTGR QVAMGLAGLV VIEDDEILKL MLPKQWGIDD VPVIVQDKKF SADGQIDYQL DVMTAAVGWF GDTLLTNGAI YPQHAAPRGW LRLRLLNGCN ARSLNFATSD NRPLYVIASD GGLLPEPVKV SELPVL MGER FEVLVEVNDN KPFDLVTLPV SQMGMAIAPF DKPHPVMRIQ PIAISASGAL PDTLSSLPAL PSLEGLTVRK LQLSMDPMLD MMGMQMLMEK YGDQAMAGMD HSQMMGHMGMH GNMNHMNHGG KDFDFHHANKI NGQAFDMNKP MFAAAKGQYE RWVISGVGDM MLHPFHIIHGT QFRILSENGK PPAAHRAGWK DTVKVEGNVS EVLVKFNHDA PKEHAYMAHC HLLHEHDTGM MLGFTV
<b>Source</b>	E.coli
<b>Target Names</b>	cueO
<b>Protein Names</b>	Recommended name: Blue copper oxidase CueO Alternative name(s): Copper efflux oxidase
<b>Expression Region</b>	29-516
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.