



Recombinant Escherichia coli Uncharacterized protein yacH (yacH)

Product Code	CSB-BP339715ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36682
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKMTLPFKPH VLALICSAGL CAASTGLYIK SRTVEAPVEP QSTQLAVSDA AAVTFPATVS APPVTPAVVK SAFSTAQIDQ WVAPVALYPD ALLSQVLMAS TYPTNVAQAV QWSHDNPLKQ GDAAIQAVSD QPWDASVKSL VAFPQLMALM GENPQWVQNL GDAFLAQPD VMDSVQRLRQ LAQQTGSLKS STEQKVITTT KKAVPVKQTV TAPVIPSNTV LTANPVITEP ATTVISIEPA NPDVVYIPNY NPTVVYGNWA NTAYPPVYLP PPAGEPFVDS FVRGFGYSMG VATTYALFSS IDWDDDDHDH HHHNDNDYHH HDGGHRDGNG WQHNGDNINI DVNNFNRTG EHLTDKNMAW RHNPNYRNGV PYHDQDMAKR FHQTDVNGGM SATQLPAPTR DSQRQAAANQ FQQRTHAAPV ITRDTQRQAA AQRFNAAEHY GSYDDFHDFS RRQPLTQQQK DAARQRYQSA SPEQRQAVRE RMQTNPKIQQ RREAARERIQ SASPEQRQAV REKMQTNPQN QQRDAARER IQSASPEQRQ VFKEKVQQR LNQQQRDNAR QRVQSASPEQ RQVFREKVQE SRPQRLNDSN HTVRLNNEQR SAVCERLSER GARRLER
Source	Baculovirus
Target Names	yacH
Protein Names	Recommended name: Uncharacterized protein yacH
Expression Region	1-617
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.