



Recombinant *Saccharomyces cerevisiae* Aminopeptidase Y (APE3)

Product Code	CSB-BP339800SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37302
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	IPKP HIPYFMKPHV ESEKLQDKIK VDDLNATAWD LYRLANYSTP DYGHPTRVIG SKGHNKTMEY ILNVFDDMQD YYDVSLQEFE ALSGKIISFN LSDAETGKSF ANTTAFALSP PVDGFGKLV EIPNLGCEEK DYASVPPRH NEKQIALIER GKCPFGDKSN LAGKFGFTAV VIYDNEPKSK EGLHGTLGEP TKHTVATVGV PYKVGKKLIA NIALNIDYSL YFAMDSYVEF IKTNIIADT KHGDPDNIVA LGAHSDSVEE GPGINDDGS G TISLLNVAKQ LTHFKINNKV RFAWWAAEEE GLLGSNFYAY NLTKEENSKI RVFMDYDMMMA SPNYEYEIYD ANNKENPKGS EELKNLYVDY YKAHHLNYTL VPFDRSDYV GFINNGIPAG GIATGAEKNN VNNGKVLDR C YHQLCDDVSN LSWDAFITNT KLIAHSVATY ADSFEGFPKR ETQKHKEVDI LNAQQPQFKY RADFLII
Source	Baculovirus
Target Names	APE3
Protein Names	Recommended name: Aminopeptidase Y EC= 3.4.11.15
Expression Region	57-537
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.