



Recombinant Escherichia coli Uncharacterized protein yhiR (yhiR)

Product Code	CSB-EP339847ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37634
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MLSYRHSFHA GNHADVLKHT VQSLIIESLK EKDKPFLYLD THAGAGRYQL GSEHAERTGE YLEGIARIWQ QDDLPAELEA YINVVKHFNR SGQLRYYPGS PLIARLLLRE QDSLQLTELH PSDYPLLRSE FQKDSRARVE KADGFQQLKA KLPPVSRRL ILIDPPYEMK TDYQAVVSGI AEGYKRFATG IYALWYPVVL RQQIKRMIHD LEATGIRKIL QIELAVLPDS DRRGMTASGM IVINPPWKLE QQMNNVLPWL HSKLVPAGTG HATVSWIVPE
Source	E.coli
Target Names	rlmJ
Protein Names	Recommended name: Uncharacterized protein yhiR
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.