



# Recombinant Escherichia coli Transcriptional regulatory protein rtcR (rtcR)

<b>Product Code</b>	CSB-EP339904ENV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P38035
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRKTVAFGFV GTVLDYAGRG SQRWSKWRPT LCLCQQESLV IDRLELLHDA RSRSLFETLK RDIASVSPET EVVSVEIELH NPWDFEEVYA CLHDFARGYE FQPEKEDYLI HITTGTHVAQ ICWFLLAEAR YLPARLIQSS PPRKKEQPRG PGEVTIIDLD LSRYNIAISR FAEERQQTL DFLKSGIATR NPHFNMIEQI EKVAIKSRAP ILLNGPTGAG KSFLARRILE LKQARHQFSG AFVEVNCATL RGDTAMSTLF GHVKGAF TGA RESREGL LRS ANGGMLFLDE IGELGADEQA MLLKAIEEKT FYPFGSDRQV S SDFQLIAGT VRDLRQLVAE GKFREDLYAR INLWTF TLPG LRQRQEDIEP NLDYEVERHA SLTGDSVRFN TEARRAWLAF ATSPQATWRG NFRELSASVT RMA T FATSGR ITLDVVEDEI NRLRYNWQES RPSALTALLG AEAENIDLFD RMQLEHVIAI CRQAKSLSAA GRQLFDVSRQ GKASVNDADR LRKYLARFGL TWEAVQDQHS SS
<b>Source</b>	E.coli
<b>Target Names</b>	rtcR
<b>Protein Names</b>	Recommended name: Transcriptional regulatory protein rtcR
<b>Expression Region</b>	1-532
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.