



Recombinant *Saccharomyces cerevisiae* Protein ECM13 (ECM13)

Product Code	CSB-BP339930SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38195
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNSTTPSIQD QYILASKVRS KLA KCVS VTT KNKDYNLRVL VGHANLLDKI TENVETHNAA TNALAGDPFS KGPENLSIEH IELSNANASK NDVGKEENAE KTIESEDYCD FYSSDEDPDA DTLSSSTDSED DDDYEDYDFE YDYSGGDYNK KIDMYFSFHT APNYQYLTH TNSHSEQTDEL AESTPRYNAL PATASTTEEE QDTETLDAVS LHSSAPIFRV LSRRVNGQED DSENESSSDV DDGSVPLTRF HSCPITA
Source	Baculovirus
Target Names	ECM13
Protein Names	Recommended name: Protein ECM13 Alternative name(s): Extracellular mutant protein 13
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.