



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YBR053C (YBR053C)

Product Code	CSB-BP339938SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38235
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSSVGDFFEEI ILHDLKPYHY VPGAIHSEGI TFKETGTLL WVDIFKGEVH KVEDIEQPES SHSFFSISRA NYGKNASIEY PPNPDELKES VGCIFPILDG ASQNEIKQVL FGSKFGIGKL DFSKSEWEYV ILYSECPELS TDRAYKLRSN DGNVSPDGKY IYVGLMSDFP FDLEPIGCLL RVDLLAHKIE LVWNCLLIPN AIHWDESDQK TMYVTDSLNF TIWKCPGGDL LKRDELIDVK NSNNQSFESP EPDGSIAWFS KDGKHSGLF ITVWSTSKVQ MFDLTNGKLL KEFILPEQTP RVSCCFVVGK DLFVTTANAE INDAVRTNTD KNGGCIYKIP NVLDGNVPLE STKRQPLH
Source	Baculovirus
Target Names	YBR053C
Protein Names	Recommended name: Uncharacterized protein YBR053C
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.