



Recombinant *Saccharomyces cerevisiae* Centractin (ARP1)

Product Code	CSB-MP339995SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38696
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDQLSDSYAL YNQPVIDNG SGIKAGFSG EERP KALEYC LVGNTKYDKV MLEGLQGDTF IGNN AQKLRG LLKLRYP IGH GVV EDWDSME LIWSYVLNEV LQLQNIG EHP LLITEAPMNP LKNREQMAQV LFETFDVSAL YVSNPAVLSL YASGR TTGCV VDCGEGYCST VPIYDGFALP ASMMRMDIGG ADITEQLQFQ LRKSAGVSLF SSSEREIVRT MKEKVCYLAK NIKKEEEKYL QGTQDLISTF KLPDGR CIEV GNDRYRAPEI LFSPQII GLG YDGLSDMCMQ SIWKVDL DLR KPLLSS ILS GGTTTLKGFG DRMLWDLEAL TKGTSKIKII APSEKYYTTW IGGSILTGLS TFQRLWTKKS DWLEDSTRVY SNLM
Source	Mammalian cell
Target Names	ARP1
Protein Names	Recommended name: Centractin Alternative name(s): Actin-like protein Actin-related protein 1
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.