



Recombinant Enterobacteria phage T4 Uncharacterized 38.9 kDa protein in nrdC-mobD intergenic region (y05G)

Product Code	CSB-MP340070EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39242
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	≥85% (SDS-PAGE)
Sequence	MKTVVKSYPFG SHLYGTSTPE SDVDFKEIFV PPARDILIGN VKEHMSKNTN NTSSKNTKDN IDHELYSLKY FFKLAADGEP VALDMLHTPP ELVVKSDLPD VWKFIQDNRS RFYTTNMKSY LGYVRKQASK YGVKGSRLAA LRDVLKVVNQ IPEQWVDYQE DGSIKQRRTK VEDIKHLRPE NEFCEWVFHN HEKTGPQTFY TVLGRKYQTT LSLIELKQSL NKLDAEYGER ARKAEANEGI DWKALSHACR GGLQLLEIYK TGDLVYPLQD APFILDVKLG KHPFKTVQEF LEDVVDQVEA ASTEASKNGM QQKVDMSFWD DFLERVYLEN HRSYYK
Source	Mammalian cell
Target Names	y05G
Protein Names	Recommended name: Uncharacterized 38.9 kDa protein in nrdC-mobD intergenic region
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.