



Recombinant Enterobacteria phage T4 Uncharacterized 36.3 kDa protein in nrdC-mobD intergenic region (y04O)

Product Code	CSB-EP340072EDZ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39254
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	≥85% (SDS-PAGE)
Sequence	MKTRKHYIDY FDSLITKHRD YQKGHREVIN NILRDFLDYI GWENHICKDT QNAYSHSLGS LLEWFKRSRL LSSVIAVNNV KKFMYPSYIE TNVSNDNVVT FNIINDVKRT YLEEWFSKDS KEKFASEFSH EFNNNVNMLF KHSRRLFCHG DDRTINVNVK DWVTAKFIPS SQNGPFELLI IVCAPHEIYK NLPYMKPCEA NKHNKTIRSL TYKLRTLLSK MDVVESFDDN TNYGLSLFET KVIKLDPN KFKPTKPNH GNDTMKEERE YLSTRLIEVE KQIEEHTKVL KDLTAKANGL RNAIEVLK
Source	E.coli
Target Names	y04O
Protein Names	Recommended name: Uncharacterized 36.3 kDa protein in nrdC-mobD intergenic region
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.