



Recombinant Escherichia coli Uncharacterized protein yjgL (yjgL)

Product Code	CSB-BP340083ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39336
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSKISDLNYS QHITLADNFK QKSEVLNTWR VGMNDFARIA GGQDNRRNIL SPGAFLEFLA KIFTLGYVDF SKRSNEAGRN MMAHIKSSSY SKDTNGNEKM KFYMNNPVGE RADSPKVIIE ISLSTITTMG TRQGHTAIIF PQPDGSTNRY EGKSFERKDE SSLHLITNKV LACYQSEANK KIARLLNNNQ ELNNLQKLNN LQKLNNLLKL NNIQGLNNPQ ELNNPQNLND SQELNNSQEL NSPQELNDPQ ELNNSQDLNN SKVSVTVSVD STITGLLKEP LNNALLAIRN EHLLLMPHVC DESISYLLGE KGILEEIDKL YALNDHGIDN DKVGNNEIND IKVNLSHILI DSLDDAKVNL TPVIDSILET FSKSPYINDV RILDWCFNKS MQYFDDTKKI KHACSVINHI NLRSDQSKIA ETLFFNLDKE PYKNSPELQG LIWNKLVVYV NEFNLSNREK TNLIQRLFDN VESIFNEVPV SILVNDIFMN DFFMKNPEMI NWYFPQLLKS YEGEKIYFDN LKYDLNDNDK ESNKEILKNQ PDNVIKEKLN NEYKLRFRMM QTILQSRVNV LPYINEQRLN KLNPPENLRI AIEHFGWKNR PITA
Source	Baculovirus
Target Names	yjgL
Protein Names	Recommended name: Uncharacterized protein yjgL
Expression Region	1-604
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.