



Recombinant *Drosophila melanogaster* Adult enhancer factor 1 (Aef1)

Product Code	CSB-MP340094DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39413
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MMHIKSLPHA HAAATAMSSN CDIVIVAAQP QTTIANNNNN ETVTQATHPA HMAAVQQQQQ QQQQQQQQH QQQQQSSGP PSVPPPTEL PLPFQMHL SG ISAEAHSAQ AAAMAAQAA AAQAAAAEQQ QPPPPTSHLT HLTTHSPTTI HSEHYLANGH SEHPGEGNAA VGVGGAVREP EKPFHCTVCD RRFRQLSTLT NHVKIHTGEK PYKCNVCDKT FRQSSTLTNH LKIHTGEKPY NCNFCPKHFR QLSTLANHVK IHTGEKPFEC VICKKQFRQS STLNNHIKIH VMDKVYVPVK IKTEEDEG
Source	Mammalian cell
Target Names	Aef1
Protein Names	Recommended name: Adult enhancer factor 1 Short name= AEF-1
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.