



Recombinant Bacillus subtilis Flagellar hook-basal body complex protein flhO (flhO)

Product Code	CSB-BP340148BRJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39752
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	>85% (SDS-PAGE)
Sequence	MLKGLYTATS AMIAQQRRT EMLSNNIANAN TSGYKADQGS MRAFPEMLLS RIESKSPAGT SRTEIGSVNT GVYMQELKPL FTQGSLKSTD QPTDIALIEN QVPMSAETNE KAALFYVPVQT ADGIRYSKSS TFSLNENNQL TINGRPILST DRQPITVDNE NFTVSENGTV TTNNGRTAGQI DVRMAEDVRN LKRDGNDLYS TADGNDLPSA AGNNQVAYSL KQGVSELSNV DVTSAYTEMT EAYRSFEANQ KVIQAYDKSM DKAANEIGKI
Source	Baculovirus
Target Names	flhO
Protein Names	Recommended name: Flagellar hook-basal body complex protein flhO
Expression Region	1-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.