



Recombinant *Drosophila melanogaster* La protein homolog (La)

Product Code	CSB-YP340302DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40796
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MAEVAETPSV EAQEEVAQPA EAQVLEAKNG DAKKDPAPAA EEAAGGFTKQ ERAIIRQVEY YFGDANLNRD KFLREQIGKN EDGWVPLSVL VTFKRLASLS TDLSEIVAAL NKSEEG LVEI SEDKLSLRRH PERPIPEHNE ERRKEIQERT AYAKGFPLDS QISELLDFAA NYDKVVNLTM RKHYDKPTKS YKFKGSIFLT FETKDQAKAF LEQEKIVYKE RELLRKWQVD YLKEKQEEYA QKNEKRKNKK EAKPEPAFEL PKNAIVVFEG APETSSREEI REAFEDIKDF EVAYIEFAKG ETKGSVRLTE ADAAEKYIAK VEEGKLFKFD EVSLSLRKAT EEEEEKEFIDK AIEFMKKRRD FTRNKGKRFN RKRHHGGNDHK HGGGKKARGD
Source	Yeast
Target Names	La
Protein Names	Recommended name: La protein homolog Alternative name(s): La autoantigen homolog La ribonucleoprotein
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.