



Recombinant Escherichia coli Protein traB (traB)

Product Code	CSB-EP340345ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41067
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MASINTIVKR KQYLWLGIVV VGTASAIGGA LYLSVDVMSG NGETVAEQEP VPDMTGVVDT TFDDKVRQHA TTEMQVTAQ MQKQYEEIRR ELDVLNKQRG DDQRRIEKLG QDNAALAEQV KALGANPVTA TGEPVPQMPA SPPGPEGEPQ PGNTPVSFPP QGSVAVPPPT AFYPGNGVTP PPQVTYQSVV VPNRIQRKVF TRNEGKQGPS LPYIPSGSFA KAMLIEGADA NASVTGNEST VPMQLRITGL VEMPNSKTYD ATGCFVGLEA WGDVSSERAI VRTRNISCLK DGKTIDMPIK GHVSFRGKNG IKGEVVMRNG KILGWAWGAG FVDGIGQGME RASQPAVGLG ATAAYGAGDV LKMGIGGGAS KAAQTLSDYY IKRAEQYHPV IPIGAGNEVT VVFQDGFQLK TVEEMALERT QSRAEEDNPE SPVPVPPSAE SHLNGFNTDQ MLKQLGNLNP QQFMSGSQGG GNDGK
Source	E.coli
Target Names	traB
Protein Names	Recommended name: Protein traB
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.