



Recombinant Escherichia coli Putative type II secretion system protein H (gspH)

Product Code	CSB-EP340386ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41443
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	FTLL EMMMLVLA LVA ITASVVLFTY GREDVA STRA RETAARFTAA LELAIDRATL SGQPVG I HFS DSAWRIMVPG KTPSAWRWVP LQEDA ADESQ NDWDEELSIH LQPFKPDDSN QPQVVILADG QITPFSLLMA NAGTGEPLLT LVCSGSWPLD QTLARDTRP
Source	E.coli
Target Names	gspH
Protein Names	Recommended name: Putative type II secretion system protein H Short name= T2SS protein H Alternative name(s): Protein transport protein HofH Putative general secretion pathway protein H
Expression Region	7-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.