



Recombinant Nuclear hormone receptor family member odr-7 (odr-7)

Product Code	CSB-BP340461CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41933
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MIVPDTEGLL IYSYGLMYGS YCMACQMLIP HFQCIPGIFP NFRISTELIK TMTDKLEQPN NNVPQQPWGP FPPAFGGRPS GEQTDGNPGE FDNDAAHQQT APFMTHFFPR IGLQFPDFTE YQRFNGFQRN AFFPNPFGSQ FTGQAFAQSF PLHNSMTTMD GFNLTHAPHP FSTNTNSTKP KDIENTVQST IKHSSENIQD KPPVLSVEYP VKYDSELKFD ANVDFTAVPK QESSDDSTLK NLKKSDQQLQ QPQQFTFPPP LLAEKSFEQP RMREDVLPFH PQFYAPLDM GTNFKQEMRT PPIDGHIDYR KFDASGKRME FQPPGALHDC QVCLSTHANG LHFGARTCAA CAAFFRRTIS DDKRYVCKRN QRCNNASRDG TGYRKICRSC RMKRCLEIGM LPENVQHKN RRDSGSPPRK TPFDTFFNGF YPSFQPSGSA AQPITVSSSE SPRHTTN
Source	Baculovirus
Target Names	odr-7
Protein Names	Recommended name: Nuclear hormone receptor family member odr-7 Alternative name(s): Odorant response abnormal protein 7
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.