



Recombinant Human 26S proteasome non-ATPase regulatory subunit 14 (PSMD14)

Product Code	CSB-YP018904HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00487
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDRLLRLGGG MPGLGQGPPT DAPAVDTAEQ VYISSLALLK MLKHGRAGVP MEVMGLMLGE FVDDYTVRVI DVFAMPQSGT GVSVEAVDPV FQAKMLDMLK QTGRPEMVVG WYHSHPGFGC WLSGVDINTQ QSFEALSERV VAVVVDPIQS VKGKVVIDAF RLINANMMVL GHEPRQTTSN LGHLNKPSIQ ALIHGLNRHY YSITINYRKN ELEQKMLLNL HKKSWMEGLT LQDYSEHCKH NESVVKEMLE LAKNYNKAVE EEDKMTPEQL AIKNVGKQDP KRHLEEHVDV LMTSNIVQCL AAMLDTVVFK
Source	Yeast
Target Names	PSMD14
Protein Names	Recommended name: 26S proteasome non-ATPase regulatory subunit 14 EC=3.4.19.- Alternative name(s): 26S proteasome regulatory subunit RPN11 26S proteasome-associated PAD1 homolog 1
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.