



Recombinant *Saccharomyces cerevisiae* Protein URE2 (URE2)

Product Code	CSB-BP340527SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23202
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNNGNQVS NLSNALRQVN IGNRNSNTTT DQSNINFEFS TGVNNNNNNN SSNNNNVQN NNSGRNGSQN NDNENNIKNT LEQHRQQQQA FSDMSHVEYS RITKFFQEQP LEGYTLFSHR SAPNGFKVAI VLSELGFHYN TIFLDFNLGE HRAPEFVSVN PNARVPALID HGMDNLSIWE SGAILLHLVN KYYKETGNPL LWSDDLADQS QINAWLFFQT SGHAPMIGQA LHFYFHSQK IASAVERYTD EVRRVYGVVE MALAERREAL VMELDTENAA AYSAGTTPMS QSRFFDYPVW LVGDKLTIAD LAFVPWNNVV DRIGINIKIE FPEVYKWTKH MMRRPAVIKA LRGE
Source	Baculovirus
Target Names	URE2
Protein Names	Recommended name: Protein URE2
Expression Region	2-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.