



Recombinant Human Nuclear factor 1 B-type (NFIB)

Product Code	CSB-BP015755HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00712
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MYSPICLTQ DEFHPFIEAL LPHVRAIAYT WFNLQARKRK YFKKHEKRMS KDEERAVKDE LLSEKPEIKQ KWASRLAKL RKDIRQEYRE DFVLTVTGKK HPCCVLSNPD QKGKIRRIDC LRQADKVVRL DLVMVILFKG IPLESTDGER LMKSPHCTNP ALCVQPHHIT VSVKELDLFL AYYVQEQDSG QSGSPSHNDP AKNPPGYLED SFVKSGVFNV SELVRVSRTP ITQGTGVNFP IGEIPSQPY HDMNSGVNLQ RSLSSPPSSK RPKTISIDEN MEPSTGDFY PSPSSPAAGS RTWHERDQDM SSPTTMKKPE KPLFSSASPQ DSSPRLSTFP QHHHPGIPGV AHSVISTRTP PPPSPLPFPT QAILPPAPSS YFSHTIRYP PHLNPQDTLK NYVPSYDPSS PQTSQSWYLG
Source	Baculovirus
Target Names	NFIB
Protein Names	Recommended name: Nuclear factor 1 B-type Short name= NF1-B Short name= Nuclear factor 1/B Alternative name(s): CCAAT-box-binding transcription factor Short name= CTF Nuclear factor 1/B Short name= NF-1/B Short
Expression Region	2-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.