



Recombinant Momordica charantia Ribonuclease MC

Product Code	CSB-BP340549MHQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23540
Product Type	Recombinant Protein
Immunogen Species	Momordica charantia (Bitter melon) (Balsam pear)
Purity	>85% (SDS-PAGE)
Sequence	FDSFWFVQQW PPAVCSFQKS GSCPGSGLRT FTIHGLWPQG SGTSLTNCPQ GSPFDITKIS HLQSQLNTLW PNVLRANNQQ FWSHEWTKHG TCSESTFNQA AYFKLAVDMR NNYDIIGALR PHAAGPNGRT KSRQAIKGFL KAKFGKFPGL RCRTDPQTKV SYLVQVVACF AQDGSTLIDC TRDTCGANFI F
Source	Baculovirus
Protein Names	Recommended name: Ribonuclease MC Short name= RNase MC EC= 3.1.27.1
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.