



Recombinant Drosophila melanogaster Protein odd-skipped (odd)

Product Code	CSB-YP340573DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23803
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSSTSASPIS NITVDDELNL SREQDFAEED FIVIKEERET SLSPMLTPPH TPTEEPLRRV HPAISEEAVA TQLHMRHMAH YQQQQQQQQQ QQQHRLWLQM QQQQQQHQP QQYPVYPTAS ADPVAVHQQQL MNHWRNAAI YQQQQQQQQH PHHHHHHGHP HHPHPHPHHV RPYPAGLHSL HAAVMGRHFG AMPTLKLGGG GGASGVPSGA TGSSRPKKQF ICKYCNRQFT KSYNLLIHER THTDERPYSC DICGKAFRRQ DHLRDHRYIH SKDKPFKCSG CGKGFCQSRT LAVHKVTHLE EGPHKCPICQ RSFNQRANLK SHLQSHSEQS TKEVVVTTSP ATSHSVPNQA LSSQPENLA QHLPLVLDLSS SSSSSEKPKR MLGFTIDEIM SR
Source	Yeast
Target Names	odd
Protein Names	Recommended name: Protein odd-skipped
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.