



Recombinant Escherichia coli CDP-diacylglycerol--serine O-phosphatidyltransferase (pssA)

Product Code	CSB-EP340575ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23830
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	<p> MLSKFKRNKH QQHLAQLPKI SQSVDDVDFF YAPADFRETL LEKIASAKQR ICIVALYLEQ DDGGKGILNA LYEAKRQRPE LDVRVLVDWH RAQRGRIGAA ASNTNADWYC RMAQENPGVD VPVYGVPIINT REALGVLHFK GFIIDDSVLY SGASLNDVYL HQHDKYRYDR YHLIRNRKMS DIMFEWVTQN IMNGRGNRL DDVNRPKSPE IKNDIRLFRQ ELRDAAYHFQ GDADNDQLSV TPLVGLGKSS LLNKTIFHLM PCAEQKLTIC TPYFNLPAIL VRNIIQLLRE GKKVEIIVGD KTANDFYIPE DEPFKIIGAL PYLYEINLRR FLSRLQYYVN TDQLVVRLWK DDDNTYHLKG MWVDDKWMLI TGNNLNPRAW RLDLENAILI HDPQLELAPQ REKELELIRE HTTIVKHYRD LQSIADYPVK VRKLIRRLRR IRIDRLISRI L </p>
Source	E.coli
Target Names	pssA
Protein Names	Recommended name: CDP-diacylglycerol--serine O-phosphatidyltransferase EC= 2.7.8.8 Alternative name(s): Phosphatidylserine synthase
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.