



Recombinant Escherichia coli Putative protein rhsE (rhsE)

Product Code	CSB-EP340604ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24211
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSTRLTTLSSH TSEGHRVSVH YGYDDKGRLT GERQTVENPE TGELLWHHET GHAYNEQGLA NRVTPDSLPP VEWLTYGSGY LAGMKLGGTP LLEFTRDR LH RETVRSFGSM AGSNAAYKLT STYTPAGQLQ SQHLNSLVYD RDYGWNDNGD LVRISGPRQT REYGYSATGR LESVRTLAPD LDIRIPYATD PAGNRLPDPE LHPDSTLTVW PDNRIAEDAH YVYRHDEYGR LTEKTDRI PA GVIRTDERT HHYHYDSQHR LVFYTRIQHG EPLVESRYLY DPLGRRMAKR VWRRELDLTG WMSLSRKPEV TWYGWGDRL TTVQDTTRI QTVYEPGSFT PLIRVETENG EREKAQRSL AETLQQEGSE NGHGVVFP AE LVRLDRLEE EIRADRV SSE SRAWLAQCGL TVEQLARQVE PEYTPARKVH FYHCDHRGLP LALISEDGNT AWRGEYDEWG NQLNEENPHH LHQPYRLPGQ QHDEESGLYY NRHRHYDPLQ GRYITPDP IGRGGWNMYQY PLNPIQVIDP MGLDAIENMT SGGLIYAVSG VPGLIAANSI TNSAYQFGYD MDAIVGGAHN GAADAMRHCY LMC RMTKTFG STIADVIGKN HEAAGDRQGG PAKERIMDLK NNTVGIACGD FSAKCSDACI EKYNTGQLFG LDGIKADNPI KAKQGSSDAS NY
Source	E.coli
Target Names	rhsE
Protein Names	Recommended name: Putative protein rhsE
Expression Region	1-682
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.