



# Recombinant Escherichia coli Putative protein rhsE (rhsE)

<b>Product Code</b>	CSB-YP340604ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P24211
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSTRLTTLSSH TSEGHRVSVH YGYDDKGRLT GERQTVENPE TGELLWHHET GHAYNEQGLA NRVTPDSLPP VEWLTYGSGY LAGMKLGGTP LLEFTRDR LH RETVRSFGSM AGSNAAYKLT STYTPAGQLQ SQHLNSLVYD RDYGWNDNGD LVRISGPRQT REYGYSATGR LESVRTLAPD LDIRIPYATD PAGNRLPDPE LHPDSTLTVW PDNRIAEDAH YVYRHDEYGR LTEKTDRI PA GVIRTDERT HHYHYDSQHR LVFYTRIQHG EPLVESRYLY DPLGRRMAKR VWRRE RDLTG WMSLSRKPEV TWYGWGDRL TTVQDTTRI QTVYEPGSFT PLIRVETENG EREKAQRSL AETLQQEGSE NGHGVVFP AE LVRLDRLEE EIRADRV SSE SRAWLAQCGL TVEQLARQVE PEYTPARKVH FYHCDHRGLP LALISEDGNT AWRGEYDEWG NQLNEENPHH LHQP YRLPGQ QHDEESGLYY NRHRHYDPLQ GRYITPDP IGRGGWNMYQY PLNPIQVIDP MGLDAIENMT SGGLIYAVSG VPGLIAANSI TNSAYQFGYD MDAIVGGAHN GAADAMRHCY LMC RMTKTFG STIADVIGKN HEAAGDRQGG PAKERIMDLK NNTVGIACGD FSAKCSDACI EKYNTGQLFG LDGIKADNPI KAKQGSSDAS NY
<b>Source</b>	Yeast
<b>Target Names</b>	rhsE
<b>Protein Names</b>	Recommended name: Putative protein rhsE
<b>Expression Region</b>	1-682
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.