



Recombinant Zea mays Zein-alpha PMS1 (ZMPMS1)

Product Code	CSB-EP340625ZAX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24449
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	TIFPQCSQA PIASLLPPYL SPAVSSVCEN PILQPYRIQQ AIAAGILPLS PLFLQQSSAL LQQLPLVHLL AQNIRAQQLQ QLVLANVAAY SQQQQFLPFN QLAALNSAAY LQQQQLLPFS QLTAAYPQQF LPFNQLAALN SAAYLQQQQL LPFSQLAVVS PAAFLTQQQL LPFYLHAVPN AGTLLQLQQL LPFNQLALTN PAAFYQQPII GGALF
Source	E.coli
Target Names	ZMPMS1
Protein Names	Recommended name: Zein-alpha PMS1 Alternative name(s): 19 kDa zein PMS1
Expression Region	22-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.