



Recombinant *Saccharomyces cerevisiae* 20S rRNA accumulation protein 4 (TSR4)

Product Code	CSB-EP340667SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25040
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSKIEELPPS DTDDHSYSSK PGDVFLAFVD APVKETDDIL VEDSFIGGEP KWLHPDSEPP AELLKCGACK SADNMKLLQ AFSPLDDEQM SAIQQRLGIN NMSYINPQDD RVLYVFLCTE CQRKGNSVRC IRGVKKNKNV DSLSEKMAST SLEKDFQINP FDLSNNSDSK CNAFSSNPFQ GANANPFGAD SINSNISQSK DEGKKKESAT VSAKTARKLH DLQKDKEYDG NCKCFKSCLLY VEEETFKNKK PAHLQLPKNL KIDKEALDLT GDEDLEKDPI KLDPRTEKLS KFLDDDDTFQK FQEVVGYNPL QVLRDLDGGK PLYAETKVD ILSTVPRPGY NPSSQRIFEM QLMPKMIFDL EEVVSVDNGM EWGTILVFTD VENYMPEFDE HGVGYVEECV KVQWESRT
Source	E.coli
Target Names	TSR4
Protein Names	Recommended name: 20S rRNA accumulation protein 4
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.