



# Recombinant Escherichia coli Arylsulfatase (asIA)

<b>Product Code</b>	CSB-EP340708ENV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P25549
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ADTPST ATARKGFAGY DHPNQYLVPK ATTIADNMMP VMQHPAQDKE TQQKLAELEK KTGKKPNVVV FLLDDVGWMD VGFNGGGVAV GNTPDIDAV ASQGLIL TSA YSQPSSSPTR ATILTGQYSI HHGILMPPMY GQPGGLQGLT TLPQLLHDQG YVTQAIGKWH MGENKESQPQ NVGFDDFRGF NSVSDMYTEW RDVHVNPEVA LSPDRSEYIK QLPFSKDDVH AVRGGEQQAI ADITPKYMED LDQRWMDYGV KFLDKMAKSD KPFFLYYGTR GCHFDNYPNA KYAGSSPART SYGDCMVEMN DVFANLYKTL EKNGQLDNTL IVFTSDNGPE AEVPPHGRTP FRGAKGSTWE GGVVRVPTFVY WKGMIQPRKS DGIVDLADLF PTALDLAGHP GAKVANLVPK TTFIDGVDQT SFFLGTNGQS NRKAEHYFLN GKLAAVRMDE FKYHVLIQQP YAYTQSGYQG GFTGTMQTA GSSVFNLYTD PQESDSIGVR HIPMGVPLQT EMHAYMEILK KYPPRAQIKS D
<b>Source</b>	E.coli
<b>Target Names</b>	asIA
<b>Protein Names</b>	Recommended name: Arylsulfatase Short name= AS EC= 3.1.6.1 Alternative name(s): Aryl-sulfate sulphohydrolase
<b>Expression Region</b>	25-551
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.