



Recombinant Human Zinc finger protein 200 (ZNF200)

Product Code	CSB-MP026574HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P98182
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMAAKVVPMP PKPKQSFILR VPPDSKLGQD LLRDATNGPK TIHQLVLEHF LTFLPKPSLV QPSQVKETL VIMKDVSSSL QNRVHPRPLV KLLPKGQVQKE QETVSLYLKA NPEELVVFED LNVFHCQEEC VSLDPTQQLT SEKEDDSSVG EMMLLAVNGS NPEGEDPERE PVENEDYREK SSDDDEMDS LVSQQPPDNQ EKERLNTSIP QKRKMRNLLV TIENDTPLEE LSKYVDISII ALTRNRRTRR WYTCPLCGKQ FNESSYLISH QRTHTGEKPY DCNHCGKSFN HKTNLNKHER IHTGEKPYSC SQCGKNFRQN SHRSRHEGIIH IREKIFKCPE CGKTFPKNEE FVLHLQSHEA ERPYGCKKCG RRFGRLSNCT RHEKTHSACK TRKQK
Source	Mammalian cell
Target Names	ZNF200
Protein Names	Recommended name: Zinc finger protein 200
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.