



Recombinant *Saccharomyces cerevisiae* 60S ribosomal protein L16-A (RPL16A)

Product Code	CSB-EP340813SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26784
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SVEPVVVID GKGHLVGRLA SVVAKQLLNG QKIVVVRAEE LNISGEFFRN KLKYHDFLRK ATAFNKTRGP FHFRAPSRIF YKALRGMVSH KTARGKAALE RLKVFEGIPP PYDKKKRVVV PQALRVLRK PGRKYTTLGK LSTSVGWKYE DVVAKLEAKR KVSSAEYYAK KRAFTKKVAS ANATAAESDV AKQLAALGY
Source	E.coli
Target Names	RPL16A
Protein Names	Recommended name: 60S ribosomal protein L16-A Alternative name(s): L13a L21 RP22 YL15
Expression Region	2-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.