



# Recombinant *Saponaria officinalis* Ribosome-inactivating protein saporin-2 (SAP2)

<b>Product Code</b>	CSB-YP340891SCG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P27559
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saponaria officinalis</i> (Common soapwort) ( <i>Lychnis saponaria</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	VTSITL DLVNPTAGQY SSFVDKIRNN VKDPNLKYGG TDIIVIGPPS KDKFLRINFQ SSRGTVSLGL KRDONLYVVAY LAMDNTNVNR AYYFKSEITS AELTALFPEA TTANQKALEY TEDYQSIEKN AQITQGDKSR KELGLGIDLL LTFMEAVNKK ARVVKNEARF LLIAIQMTAE VARFRYIQNL VTKNFPNKFD SDNKVIQFEV SWRKISTAIY GDAKNGVFNK DYDFGFGKVR QVKDLQMGLL MYLGPKKSSN EANSTAYATT VL
<b>Source</b>	Yeast
<b>Target Names</b>	SAP2
<b>Protein Names</b>	Recommended name: Ribosome-inactivating protein saporin-2 Short name= SAP-2 Short name= SO-2 EC= 3.2.2.22 Alternative name(s): rRNA N-glycosidase
<b>Expression Region</b>	25-292
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.