



Recombinant Escherichia coli UDP-N-acetyl-D-mannosamine dehydrogenase (wecC)

Product Code	CSB-YP340917ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27829
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSFATISVIG LGYIGLPTAA AFASRQKQVI GVDINQHAVD TINRGEIHIV EPDLASVVKT AVEGGFLRAS TTPVEADAWL IAVPTPFKGD HEPDMTYVES AARSIAPVLK KGALVILEST SPVGSTEKMA EWLAEMRPDL TFPQQVGEQA DVNIAYCPER VLPQQVMVEL IKNDRVIGGM TPVCSARASE LYKIFLEGEC VVTNSRTAEM CKLTENSFRD VNIAFANELS LICADQGINV WELIRLANRH PRVNILQPGP GVGGHCIAVD PWFIVAQNPQ QARLIRTARE VNDHKPFWVI DQVKAADVDC LAATDKRASE LKIACFGLAF KPNIDDLRES PAMEIAELIA QWHSGETLVV EPNIHQLPKK LTGLCTLAQL DEALATADVL VMLVDHSQFK VINGDNVHQQ YVVDAGVWR
Source	Yeast
Target Names	wecC
Protein Names	Recommended name: UDP-N-acetyl-D-mannosamine dehydrogenase EC=1.1.1.n3 Alternative name(s): UDP-ManNAc 6-dehydrogenase
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.