



Recombinant Arabidopsis thaliana Trichome differentiation protein GL1 (GL1)

Product Code	CSB-EP340923DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27900
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRIRRRDEKE NQEYKKGLWT VEEDNILMDY VLNHGTGQWN RIVRKTGLKR CGKSCRLRWM NYLSPNVNKG NFTEQEEDLI IRLHKLLGNR WSLIAKRVP RTDNQVKNYW NTHLSKCLVG DYSSAVKTTG EDDDSPPSLF ITAATPSSCH HQQENIYENI AKSFNGVVSA SYEDKPKQEL AQKDVLMATT NDP SHYYGNN ALWVHDDDFE LSSLVMMNFA SGDVEYCL
Source	E.coli
Target Names	GL1
Protein Names	Recommended name: Trichome differentiation protein GL1 Alternative name(s): Myb-related protein 0 Short name= AtMYB0 Protein GLABRA 1 Protein GLABROUS 1 R2R3-MYB transcription factor GL1
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.