



Recombinant Abrus precatorius Abrin-c, Partial

Product Code	CSB-MP340994AAC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28590
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Abrus precatorius (Indian licorice) (Glycine abrus)
Purity	>85% (SDS-PAGE)
Sequence	QDQVIK FTTEGATSQS YKQFIEALRQ RLTGGLIHDI PVLDPPTTVE ERNRYITVEL SNSERESIEV GIDVTNAYVV AYRAGSQSYF LRDAPASAST YLFPGTQRYS LRFDGSGYDL ERWAHQTREE ISLGLQALTH AISFLRSGAS NDEEKARTLI VIIQMASEAA RYRYISNRVG VSIRTGTAFQ PDPAMLSLEN NWDNLSGGVQ QSVQDTPNN VILSSINRQP VVVDLSHPT VAVLALMLFV CNPPN
Source	Mammalian cell
Target Names	N/A
Protein Names	Recommended name: Abrin-c Cleaved into the following 3 chains: 1. Abrin-c A chain EC= 2. 3.2.2.22 Alternative name(s): rRNA N-glycosidase Linker peptide Abrin-c B chain
Expression Region	35-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.