



# Recombinant FBR murine osteosarcoma virus Gag polyprotein (gag)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP341041FAD   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P29175  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | FBR murine osteosarcoma virus (FBR-MSV) (Finkel-Biskis-Reilly murine osteosarcoma virus)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | GQVTTPPLS LTLEHWGDVQ RIASNQSV DV KKRRWVTFCS AEWPTFDVGV<br>PQDGTFNLDI ILQVKSKVFS PGPHGHPDQV PYIVTWEAIA YEPPPWVKPF<br>VSPKLSPTSPT APILPSGPST QPPRSALYP ALTPSIKPRP SKPQVLSDDG<br>GPLIDLLTED PPPYGEQGPS SSDGDGDREE ATSTSEIPAP SPMVSRLRGK<br>RDPPAADSTT SRAFPLRLGG NGQKNNNPSF SEDPGKLTAL IESVLTTTHQP<br>TWDDCQQLLG TLLTGEEKQR VLLEARKAVR GNDGRPTQMP NEVNAAFPLE<br>RPDWDYTTPE                                     |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | gag   |
| <b>Protein Names</b>     | Recommended name: Gag polyprotein Alternative name(s): Core polyprotein<br>Cleaved into the following 3 chains: 1. Matrix protein p15 Short name= 2. MA 3.<br>RNA-binding phosphoprotein p12 Alternative name(s): pp12 Capsid protein   |
| <b>Expression Region</b> | 2-310   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |