



Recombinant Alkanal monooxygenase beta chain (luxB)

Product Code	CSB-MP341048EUW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29239
Product Type	Recombinant Protein
Immunogen Species	Photobacterium leiognathi
Purity	>85% (SDS-PAGE)
Sequence	MNFGFLFFLN F QPEGMTSEMV LDNMVDTVAL VDKDDYHFKR VLVSEHHFSK NGIIGEPLTA ISFLLGLTKR IEIGSLNQVI TTHHPVRIGE QTGLLDQMSY GRFVLGLSDC VNDFEMDFFK RKRSSQQQQF EACYEILNEA LTTNYCQADD DFFNFPRISV NPHCISEVKQ YILASSMGVV EWAARKGLPL TYRWSDSLAE KEKYYQRYLA VAKENNIDVS NIDHQFPLL V NINENRRIAR DEVREYIQSY VSEAYPTDPN IELRVEELIE QHAVGKVDEY YDSTMHAVKV TGSKNLLLSF ESMKNKDDVT KLINMFNQKI KDNLIK
Source	Mammalian cell
Target Names	luxB
Protein Names	Recommended name: Alkanal monooxygenase beta chain EC= 1.14.14.3 Alternative name(s): Bacterial luciferase beta chain
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.