



Recombinant Solanum tuberosum Actin-101 (AC101)

Product Code	CSB-EP341127FIG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30173
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MADAEDIEPL VCDNGTGMVK AGFAGDDAPR AVFPSIVGRP RHTGVMVGMG QKDAYVGDEA QSKRGILTLK YPIEHGIVSN WDDMEKIWHH TFYNELRVSP DEHPVLLTEA PLNPKANREK MTQIMFETFN VPAMYVAIQA VLSLYASGRT TGIVLDSGDG VSHTVPIYEG YALPHAILRL DLAGRDLTDC LMKILTERGY SFTTSAEREI VRDMKEKLAY VALDYEQELE TAKSSSAVEK SYELPDGQVI TIGAERFRCP EVLFQPSLVG MEAAGIHETT YNSIMKCDVD IRKDLYGNIV LSGGTTMFPG IADRMSKEIT ALAPSSMKIK VVAPPERKYS VWIGGSILAS LSTFQQMWIT KGEYDESGPS IVHRKCF
Source	E.coli
Target Names	AC101
Protein Names	Recommended name: Actin-101
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.