



Recombinant Escherichia coli Putative molybdenum transport ATP-binding protein ModF (modF)

Product Code	CSB-BP341197ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31060
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MSSLQILQGT FRLSDTKTLQ LPQLTLNAGD SWAFVGSNGS GKSALARALA GELPLLKGER QSQFSHITRL SFEQLQKLVS DEWQRNNTDM LGPGEDDTGR TTAETIIQDEV KDAPRCMQLA QQFGITALLD RRFKYLSTGE TRKTLLCQAL MSEPDLILLD EFPDGLDVAS RQQLAERLAS LHQSGITLVL VLNRFDEIPE FVQFAGVLAD CTLAETGAKE ELLQQALVAQ LAHSEQLEGV QLPEPDEPSA RHALPANEPRI IVLNNGVVSY NDRPILNLS WQVNPGEHWQ IVGPNAGGKS TLLSLVTGDH PQGYSNDLTL FGRRRGSGET IWDIKKHIGY VSSSLHLDYR VSTTVRNIL SGYFDSIGIY QAVSDRQQKL VQQWLDILGI DKRTADAPFH SLSWGQQRLA LIVRALVKHP TLLILDEPLQ GLDPLNRQLI RRFVDVLISE GETQLLFVSH HAEDAPACIT HRLEFVPDGG LYRYVLTKIY
Source	Baculovirus
Target Names	modF
Protein Names	Recommended name: Putative molybdenum transport ATP-binding protein ModF Alternative name(s): Photorepair protein PhrA
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.