



Recombinant Human Glycerol kinase (GK), partial

| | |
|--------------------------|---|
| Product Code | CSB-MP341295HU |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P32189 |
| Storage Buffer | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Product Type | Recombinant Proteins |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | PLVGAVDQGTSSSTRFLVFNSKTAELLSHHQVEIKQEFPREGWVEQDPKEILHS VYECIEKTCEKLGQLNIDISNIKAIGVSNQRETTVVWWDKITGEPLYNAVWDLDR TQSTVESLSKRIPGNNNFVKSKTGLPLSTYFSAVKLRWLLDNVRKVQKAVEEK RALFGTIDSWLIWSLTGGVNGGVHCTDVTNASRTMLFNIHSLEWDKQLCEFFG IPMEILPNVRSSSEIYGLMKISHSVKAGALEGVPISGCLGDQSAALVGQMCQFI GQAKNTYGTGCFLLCNTGHKCVFSDHGLLTTVAYKLGRDKPVYYALEGSVAIA GAVIRWLRDNLGIIKTSEEIEKLAKEVGTSYGCYFVPAFSGLYAPYWEPSARGII CGLTQFTNKCHIAFAALEAVCFQTREILDAMNRDCGIPLSHLQVDGGMSTNKIL MQLQADILYIPVVKPSMPETTALGAAMAAGAAEGVGVWSLEPEDLSAVTMERF EPQINAESEIRYSTWKKAVMKSMGW |
| Source | Mammalian cell |
| Target Names | GK |
| Protein Names | Recommended name: Glycerol kinase Short name= GK Short name= Glycerokinase EC= 2.7.1.30Alternative name(s): ATP:glycerol 3- phosphotransferase |
| Expression Region | 11-519aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.