



# Recombinant *Saccharomyces cerevisiae* Protein MSS51, mitochondrial (MSS51)

<b>Product Code</b>	CSB-YP341308SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P32335
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	VVSH QTRRHLMGFV RNALGLDPPP SPEDPTPENR FHPWDQSPSV DLRERAAKIR TLAHCPVTGK DINYTCPLSG IPTHHSREAW EMDKAYHDSK KYEILKKVNI YEHDLRSGRP FPEFDFPQQQ GYDKAVNLTN WDLFFYTRSF YSMDTEFQLA AVTKMLSYPY TIGSLLHKFS PYSLNPKGPI TLEGLKSLAA LRYTLYPLEN RSLPTTTKNR AMRIFILGAR AEAQLPGHVW KQLQFLFPEQ SFEIHFIGPE CLYKRDKQEY VKSTTPVVQR VDETLKFIYR TNFFEVFHEA QDFFPYDPYM DVFFTFHPGY ASPESHGSWM GETMKALLET KCAIFTTGFN KKDLTDDINL VKSKYGKEMD VLMEPVRNVF GSTKWELNDM NPQEVYQFNM YIAGFRGKRY HTIKRQ
<b>Source</b>	Yeast
<b>Target Names</b>	MSS51
<b>Protein Names</b>	Recommended name: Protein MSS51, mitochondrial Alternative name(s): Mitochondrial splicing suppressor protein 51
<b>Expression Region</b>	27-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.