



Recombinant Bacillus cereus MreB-like protein (mbl)

Product Code	CSB-MP341321BQN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32444
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ATCC 10987 / NRS 248)
Purity	>85% (SDS-PAGE)
Sequence	MFARDIGIDL GTANVLIHVK GKGIVLNEPS VVAIDRNTGK VLAVGEEARS MVGRTPGNIV AIRPLKDGVI ADFEITEAML KYFINKLDVK SFFSKPRILI CCPTNITSVE QKAIREAAER SGGKTVFLEE EPKVAAVGAG MEIFQPSGNM VVDIGGGTTD IAVLSMGDIV TSSSIKMAGD KFDMEILNYI KRKYKLLIGE RTSEDIKIKV GTVFPGARSE ELEIRGRDMV TGLPRTITVC SEEITEALKE NAAVIVQAAK GVLERTPEL SADIIDRGVI LTGGGALLHG IDMLLAEELK VPVLIAENPM HCVAVGTGIM LENIDRLPKR ALR
Source	Mammalian cell
Target Names	mbl
Protein Names	Recommended name: MreB-like protein Short name= Protein mbl
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.